Instructions for use

1 step vitrectomy cannula system sterile

H.P. Braem AG
Industriestrasse 4
9552 Bronschhofen
Schweiz

T: +41 (0)71 866 19 33
F: +41 (0)71 866 19 80
E: info@hpbraem.ch
W: www.hpbraem.ch

Precision is our passion

Symboles

Manufacturer

Article
number

Lot
number

Expiry
date

EO
sterilized

Do not use
if damaged

Single use
only

Store away
from sunlight

Keep
dry





Attention

CE Mark
Notified Body

Application

DE-Gebrauchsanweisung

REF15550, REF15560, REF15650, REF15750

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1 step vitrectomy cannula system sterile

Intended use: Trocarsystem for vitrectomy

Contraindications: There are no specific contraindications for the VCS.

Risks: Hypotony / Wound leakage, potentially requiring intraoperative suturing / Minor subconjunctival bleeding / Vitreous incarceration / Vitreous bleeding / Subconjunctival silicone oil / Cataract / Bacterial contamination of the vitreous cavity / potentially leading to endophthalmitis / Retinal or choroidal detachment

Set contains: 3 inserters (trocars) / 3 instrument cannulas with bulb where indicated / 3 closure valves / 1 infusion cannula with infusion tubing and luer / 3 stoppers in a box (only by 23G)

Caption: A: Back of handle: Distance mark / B: Insertion angle at 20° - 30° / C: Removing the inserter / D: Infusion tubing, free removable / E: If finished: remove cannula with valve

Application: The 3 universal cannulas are to be used preloaded with the sharp stiletto inserter and closure valve. With the back of handle mark 3.5mm distance. Perform a transconjunctival/scleral incision with the preloaded sharp stiletto inserter at an angle of about 20 - 30 degrees till the middle of the cannula shaft. To overcome the somewhat greater resistance at introduction of the cannula, the inserter must be positioned perpendicular to the sclera to ensure enough counter pressure to prevent the eye from rotating away when pressure is applied. Fixate the cannula with a forceps and retract the sharp stiletto inserter. **27G Infusion** - The colored sealing valve must be removed for infusion. Stick the infusion tubing end in the first inserted cannula and switch constant infusion on. Insert the 2 other instrument cannulas in the same way as described above. The vitrectomy system is now ready to start surgery. When surgery is completed, the instrument cannulas can be removed by a forceps and finally the cannula with infusion cannula can be removed. In case of small subconjunctival haemorrhage from episcleral vessels, use a cotton wool applicator to press it gently to the area of the scleral tunnel.

May only be used by a specialist! Patient selection and choice of surgical methods are the sole responsibility of the ophthalmic surgeon.

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