Instructions for use

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EN-Instruction for use

REF17100

Silicone oil aspiration adapter sterile

Intended use: Aspiration adapter for the removal of silicone oil during vitreoretinal surgery. Suitable for use with sili-cone oils up to **5000 cSt**. For use with H. P. BRAEM AG vitrectomy cannula systems only. Not suitable for use with silicone oils heavier than-water!

Contraindications: There are no specific contraindications.

Risks: Hypotony / Wound leakage, potentially requiring intraoperative suturing / Minor subconjunctival bleeding / Vitreous incarceration / Vitreous / Subconjunctival silicone oil / Cataract / Bacterial contamination of the vitreous cavity potentially leading to endophthalmitis / Retinal or choroidal detachment / Lowered intraocular pressure (IOP) due to too fast aspiration

Set contains: 5 pcs. Silicone oil aspiration adapter, sterile

Caption: **A:** Insertion of infusion cannula: Once the 1 or 2 step vitrectomy cannula system 23G - 27G by H. P. BRAEM AG has been assembled and attached according to instructions, start by attaching the infusion cannula with infusion tubing to ensure stability of intraocular pressure / **B:** Attaching the oil aspiration system: The silicone oil aspiration adaptor can now be attached. To do this, remove the closure valve from the instrument cannula / **C:** Attaching the aspiration adapter to the cannula. Once this has been attached, use the luer lock to fit the aspiration adapter to the syringe connected to the active aspiration device.

Application: **Oil aspiration:** Once the vitrectomy system's active aspiration has been activated, the vacuum produced will allow silicone oil aspiration at up to **400mmHg / 0.53bar**. Ensure that infusion pressure remains adequate in relation to the degree of vacuum used for aspiration. **Disconnecting the oil aspiration system:** Once all the silicone oil has been removed from the eye, remove the aspiration tubing and use the closure valve to re-seal the instrument cannula. You may now proceed to the next stage of the surgical procedure. For information on how to remove the cannulas from the eye, please refer to the instruction manual of the vitrectomy cannula system used in the procedure.

Attention! May only be used by a specialist! Patient selection and choice of surgical methods are the sole responsibility of the ophthalmic surgeon.

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